

**色調段階に対応する濃度**  
(クリニテックアトラス XL : クリニテック 200+フォーマット)

項目 判定段階	ブドウ糖		ケトン体	蛋白質			潜血	ビリルビン	白血球
	JCCLS 対応	慣用単位		試薬カートリッジ -10	試薬カートリッジ PRO12				
			低濃度域		高濃度域				
-	0	0	0	0	0	---	0	0	0
+/-	---	100 (0.1)	5	15	15(Low)	---	0.015	---	15
1+	100	250 (0.25)	15	30	---	30	0.062	0.8	70
2+	250	500 (0.5)	40	100	---	100	0.135	1.6	125
3+	500	1000 (1.0)	80	300	---	300	0.405	3.2	500
4+	1000	---	---	---	---	---	---	---	---
感度 単位	75-125 mg/dL (g/dL)		5-10 mg/dL	15-30 mg/dL	8-15 mg/dL	30-65 mg/dL	0.015-0.062 mg/dL	0.4-0.8 mg/dL	5-15 個/hpf cell/ $\mu$ L

ウロビリノーゲン	
判定段階	
+/-	0.1
	1.0
1+	2.0
2+	4.0
3+	8.0
Ehrlich 単位/dL	

PH
判定段階
5.0
5.5
6.0
6.5
7.0
7.5
8.0
8.5

色調
判定段階
YELLOW/STRAW
DK YELLOW/AMBER
RED
ORANGE
GREEN
OTHER

亜硝酸塩	
判定段階	
-	0
+	0.06-0.1
mg/dL	

クレアチニン
判定段階
10
50
100
200
300
mg/dL