

US3500色調段階に対する濃度

項目 判定段階	ブドウ糖	ケトン体	タンパク質	潜血		ビリルビン	白血球	ウレリノゲン	亜硝酸塩
				ヘモグロビン	赤血球				
—	0	0	0	0	0	0	0	normal	0
+/-	50	...	15	0.03	10
1+	100	10	30	0.06	20	0.5	25	2.0	0.05~
2+	250	30	100	0.15	50	1	75	4.0	...
3+	500	80	300	0.75	250	2	500	8.0	...
4+	1000	...	1000	12.0	...
感度	40~50	8~10	12~15	0.024~0.03	8~10	0.4~0.5	20~25	1.6~2.0	0.05
単位	mg/dL	mg/dL	mg/dL	mg/dL	c/ μ L	mg/dL	c/ μ L	mg/dL	mg/dL

尿色
判定段階
LIGHT YELLOW
STRAW
YELLOW
AMBER
DARK BROWN
RED
OTHER

混濁
判定段階
—
1+
2+

pH
判定段階
5.0
5.5
6.0
6.5
7.0
7.5
8.0
8.5
9.0

クレアチニン
判定段階
10
50
100
200
300
mg/dL

アルブミン
判定段階
10
30
80
150
over
mg/L